


Accessories and packing material




SP Medical has developed a very flexible programme of packing materials to supply optimal protection of Guide Wires and Catheters.

Clipless Dispensers

The new design – No Clip Penetration – ensures the highest level of sterility during transportation and storage.

- Suitable for Guide Wires from 40 – 200 cm
- Internal Dispenser diameter from 9 to 19 cm
- Standard tubing dimension is Ø 4,0 x 2,75 mm.
- Tubing material is PP



Dispenser with clips

- Suitable for Guide Wires of any length up to 450 cm
- Internal Dispenser diameter is 14,5 cm
- Standard tubing dimension is Ø 4,0 x 2,75 mm
- Tubing material is PP or PE

Customised products on request

- Variable Dispenser shape
- Variable Dispenser diameter
- Variable tubing dimensions
- Variable tubing material and colour

Available packing components

- J-Straightner in different colours
- Female Luerlock with protection grid
- 2, 3 and 4 cavity clips
- Peel pouches
- Tip protector