



Accoat Standard Guide Wire – Lunderquist PTC



Dimensions

Ø 0.025" / 0,64 mm

Ø 0.035" / 0,89 mm

Length

- 50 cm
- 60 cm
- 70 cm
- 80 cm
- 90 cm
- 100 cm
- 120 cm
- 150 cm
- 260 cm

Surface

- Stainless steel
- PTFE coated

Distal tip form

- Straight tip
- J-tip radius: 1,5 – 3 – 6 mm

Distal flexible length

- 5 cm
- 7 cm
- 10 cm
- 12 cm
- 15 cm
- 17 cm

Customised products on request

Description of item code

PTC-S/Jx-C-Fx-0xx-L

(example: PTC-S-C-F5-035-150)

S	Straight tip
Jx	J-tip radius J1,5 mm tip J3 mm tip J6 mm tip
C	PTFE coated
Without C	stainless steel
Fx	Distal Flexible length F5 - 5 cm flexible tip F7 - 7 cm flexible tip F10 - 10 cm flexible tip F12 - 12 cm flexible tip F15 - 15 cm flexible tip F17 - 17 cm flexible tip
0xx:	Diameter of the guide wire 025: 0.025" (0,64 mm) 035: 0.035" (0,89 mm)
L:	Length of the guide wire 50, 60, 70, 80, 90, 100, 120, 150 and 260 cm