



**CHAMPION™ PTCA Guide Wire – For Cardiology use (CAR)**



**Dimension**  
Ø 0.014" / 0,36 mm

**Length**  
175 cm  
195 cm  
300 cm

**Distal tip form**  
Straight  
J-tip

**Stent zone stiffness**  
Standard  
Stiff

**Tip stiffness**  
Standard  
Intermediate  
Floppy  
Extra Floppy

**Radiopaque length**  
3 cm

**Coil surface**  
Uncoated Stainless Steel  
Hydrophilic coated

**Shaft surface**  
PTFE coated

**Description of item code**

PTCA-CAR-S/J-M/HY/SI-EF/F/I/ST-D/H-014-L  
(example: PTCA-CAR-S-HY-F-D-014-175)

<b>CAR</b>	For cardiology use
<b>S</b>	Straight tip
<b>J</b>	J-tip
<b>M</b>	Coil surface uncoated – stainless steel
<b>HY</b>	Coil surface hydrophilic
<b>SI</b>	Coil surface silicone
<b>EF</b>	Tip stiffness extra floppy
<b>F</b>	Tip stiffness floppy
<b>I</b>	Tip stiffness intermediate
<b>ST</b>	Tip stiffness standard
<b>D</b>	Stent zone stiffness standard
<b>H</b>	Stent zone stiffness stiff
<b>014:</b>	Diameter of the guide wire 0.014" (0,36 mm)
<b>L:</b>	Length of the guide wire Available in 175, 195, 300 cm