



**CHAMPION™ PTCA Guide Wire – Stiff tip for highly stenosed lesions (CTO)**



**Dimension**

Ø 0.014" / 0,36 mm

**Length**

195 cm  
300 cm

**Distal tip form**

Straight  
J-tip

**Stent zone stiffness**

Standard

**Tip stiffness (Gram Force)**

3,0gf  
4,5gf  
6,0gf  
9,0gf  
12,0gf

**Radiopaque length**

3 cm

**Coil surface**

Hydrophilic coated  
Silicone coated

**Shaft surface**

PTFE coated

**Description of item code**

PTCA-CTO-S/J-SI/HY-3/4,5/6/9/12-D-014-L  
(example: PTCA-CTO-S-HY-3-D-014-175)

**CTO** For highly stenosed lesions

**S** Straight tip

**J** J-tip

**HY** Coil surface hydrophilic

**SI** Coil surface silicone

**3** Tip stiffness 3 Gram Force

**4,5** Tip stiffness 4,5 Gram Force

**6** Tip stiffness 6 Gram Force

**9** Tip stiffness 9 Gram Force

**12** Tip stiffness 12 Gram Force

**D** Stent zone stiffness standard

**014:** Diameter of the guide wire 0.014" (0,36 mm)

**L:** Length of the guide wire Available in 195, 300 cm