



Accoat Standard Guide Wire – Schüller PTC



Dimensions

Ø 0.025" / 0,64 mm
Ø 0.035" / 0,89 mm

Length

50 cm
60 cm
70 cm
80 cm
90 cm
100 cm
110 cm
120 cm
150 cm
260 cm

Surface

Stainless steel
PTFE coated

Distal tip form

Straight tip
J-tip radius:
1,5 – 3 – 6 mm

Customised products on request

Distal flexible length

5 cm
7 cm
10 cm
12 cm
15 cm
17 cm

Proximal flexible length

30 cm
33 cm
35 cm
38 cm
40 cm
43 cm
45 cm
50 cm
53 cm
55 cm
60 cm
63 cm
65 cm

Description of item code

PTC-SCH-S/Jx-C-Fx-yy-zz-0xx-L
(example: PTC-SCH-S-F5-40-55-035-100)

S	Straight tip
Jx	J-tip radius J1,5 mm tip J3 mm tip J6 mm tip
C	PTFE coated
Without C	stainless steel
Fx	Distal Flexible length x = cm ex. F5 = 5 cm Available in 5, 7, 10, 12, 15 or 17 cm
yy	Shaft length - cm
zz	Proximal Flexible length - cm Available in 30, 33, 35, 38, 40, 43, 45, 50, 53, 55, 60, 63 or 65 cm
0xx:	Diameter of the guide wire 025: 0.025" (0,64 mm) 035: 0.035" (0,89 mm)
L:	Length of the guide wire (Sum of Fx+yy+zz) Available in 50, 60, 70, 80, 90, 100, 110, 120, 150 or 260 cm