


Accoat Guide Wire – Seldinger



Dimensions

- Ø 0.018" / 0,46 mm
- Ø 0.021" / 0,53 mm
- Ø 0.025" / 0,64 mm
- Ø 0.028" / 0,71 mm
- Ø 0.032" / 0,81 mm
- Ø 0.035" / 0,89 mm
- Ø 0.038" / 0,97 mm

Length

30 – 450 cm

Surface

Stainless steel
PTFE coated

Distal tip form

Straight tip
J-tip radius: 1,5 – 3 – 6 – 15 mm

Distal flexible length

Standard 4,5 cm flexible tip

Shaft stiffness

Standard
Stiff
Extra stiff
Super stiff
(Extra stiff & Super stiff not available in dim: 0.018"; 0.021"; 0.025")

Markers

Uncoated stainless steel
OEM
Length 30-150 cm
1 mm wide surface markers
10 cm step

Customised products on request

Description of item code

(H/X/X)-S/Jx-(C)-0xx-L
(example: J3-C-035-150)

H	Stiff Guide Wire
H-X	Extra Stiff Guide Wire
H-X-X	Super Stiff Guide Wire
S	Straight tip
Jx	J-tip radius J1.5 mm J3 mm J6 mm J15 mm
C	PTFE coated
Without C	stainless steel
M	Surface marker
0xx:	Diameter of the guide wire 018: 0.018" (0,46 mm) 021: 0.021" (0,53 mm) 025: 0.025" (0,64 mm) 028: 0.028" (0,71 mm) 032: 0.032" (0,81 mm) 035: 0.035" (0,89 mm) 038: 0.038" (0,97 mm)
L:	Length of the guide wire (from 30 – 450 cm)