


## Packing Material



SP Medical has developed a very flexible programme of packing materials to supply optimal protection of guide wires and catheters.

### Clipless Dispensers

The new design – No Clip Penetration – ensures the highest level of sterility during transportation and storage.

Suitable for guide wires from 40 – 200 cm  
Internal Dispenser diameter from 9 to 19 cm  
Standard tubing dimension is  $\varnothing$  4.0 x 2.75 mm.  
Tubing material is PP

### Dispenser with clips

Suitable for guide wires of any length up to 450 cm  
Internal Dispenser diameter is 14.5 cm  
Standard tubing dimension is  $\varnothing$  4.0 x 2.75 mm  
Tubing material is PP or PE

### Customised products on request

Variable Dispenser shape  
Variable Dispenser diameter  
Variable tubing dimensions  
Variable tubing material and colour

### Available packing components

J-Straightner in different colours  
Female Luerlock with protection grid  
2, 3 and 4 cavity clips  
Peel pouches  
Tip protector