



Accoat Guide Wire – Seldinger



**Dimensions**

- Ø 0.018" / 0.46 mm
- Ø 0.021" / 0.54 mm
- Ø 0.025" / 0.64 mm
- Ø 0.028" / 0.72 mm
- Ø 0.032" / 0.82 mm
- Ø 0.035" / 0.89 mm
- Ø 0.038" / 0.97 mm

**Length**

20 – 450 cm

**Surface**

- Stainless steel
- PTFE coated

**Distal tip form**

- Straight tip
- J-tip radius: 1,5 – 3 – 6 – 15 mm

**Distal flexible length**

Standard 4,5 cm flexible tip

**Shaft stiffness**

- Standard
  - Stiff
  - Extra stiff
  - Super stiff
- (Extra stiff & Super stiff not available in dim: 0.018"; 0.021"; 0.025")*

**Markers**

- Uncoated stainless steel
- OEM
- Length 30-150 cm
- 1 mm wide surface markers
- 10 cm step

**Customised products on request**

**Description of item code**

(H/X/X)-S/Jr-(C)-0xx-L  
*(example: J3-C-035-150)*

<b>H</b>	Stiff guide wire
<b>H-X</b>	Extra Stiff guide wire
<b>H-X-X</b>	Super Stiff guide wire
<b>S</b>	Straight tip
<b>Jr</b>	J-tip radius J1.5 mm J3 mm J6 mm J15 mm
<b>C</b>	PTFE coated
<b>Without C</b>	Stainless steel
<b>M</b>	Surface marker
<b>0xx</b>	Diameter of the guide wire <b>018:</b> 0.018" (0.46 mm) <b>021:</b> 0.021" (0.54 mm) <b>025:</b> 0.025" (0.64 mm) <b>028:</b> 0.028" (0.72 mm) <b>032:</b> 0.032" (0.82 mm) <b>035:</b> 0.035" (0.89 mm) <b>038:</b> 0.038" (0.97 mm)
<b>L</b>	Length of the guide wire ■ 20 – 450 cm